



**AUG**  
MEDICAL

Quality. Defined

# VETERINARY CATALOGUE





**AUG**  
MEDICAL

Quality. Defined

# WHY US



AUG Medical is an accredited manufacturer of high-quality Laryngoscopes, Diagnostics, Airway Management, Video Laryngoscopes and Otoscopes.

We offer an extensive range of Veterinary ENT devices for Animals. These devices are crafted with High Grade Medical Stainless Steel. Our outstanding wide range of Animal ENT Devices are remarkably durable and reliable as per the needs of veterinarians.

We are the direct supplier and distributor of medical equipments based in Poland, United States, UAE & Turkey. Our headquarter is in Poland which serve the whole Europe. US captures North and South America, and Turkey & UAE controls Asian Pacific. Our goal is to provide Our Customer certified and standardized products with premium quality at competitive prices with maintained stocks we provide fast and efficient delivery for the best satisfaction of our customers.

Convenient customer service is very important to us, that's why our customer service is available in more than four languages to answer to your questions or redress any concern you may have. Moreover to make purchasing more convenient for our valuable customers, we offer online shopping portal for both end-users & distributors. Our authorised distributors are listed on our website for your ease. Shopping in our online store is easy and safe because of multipurpose payment facilities & best delivery services. We also constantly improve our offers in order to satisfy our customers and provide access to the best product. The confidence offer to our brand and the full satisfaction of our customers is the answer to effective quality management and internal audits of ISO 9001 and ISO 13485. We are also accredited with cGMP, CE mark, REACH and FDA.

Thank you for choosing AUG MEDICAL.



ISO 9001  
ISO 13485



**AUG**  
MEDICAL

Quality. Defined



# INDEX

## 1. VETSKOP®

### LARYNGOSCOPE

PG No.

VETSKOP®   F.O RABBIT BLADE	02
VETSKOP®   F.O LEFT HAND MACINTOSH BLADE	02
VETSKOP®   F.O LEFT HAND FLEXI-TIP BLADE	03
VETSKOP®   F.O LEFT HAND MILLER BLADE	03
VETSKOP®   F.O LEFT HAND MILLER BLADE WITH MAGNIFICATION	04
VETSKOP®   F.O LEFT HAND MILLER BLADE(BAYONET FITTING) WITH MAGNIFICATION	04
VETSKOP®   F.O LED/XENON HANDLES	05
VETSKOP®   POWER SOURCE RHEOSTAT HANDLE(BAYONET)	05

## 1. VETSKOP®

### DIAGNOSTICS

VETSKOP®   F.O LED OTOSCOPE	07
VETSKOP®   F.O XENON OTOSCOPE	07
VETSKOP®   F.O VACUUM OTOSCOPE	08
VETSKOP®   F.O LED OPHTHALMOSCOPE	08
VETSKOP®   F.O XENON OPHTHALMOSCOPE	09
VETSKOP®   F.O RC LED/XENON OPHTHALMOSCOPE	09
VETSKOP®   OTO-OPHTHALMOSCOPE SETS & DIAGNOSTICS	10
VETSKOP®   ACCESSORIES	11

# VETSKOP® Series

## LARYNGOSCOPE



## DELINEATION OF VETSKOP® SERIES

AUG VETSKOP® Series are designed with enclosed fiber optic bundle with a whole new experience. These blades has no cavities to trap dirt or body fluids thus allowing the blade to be easily cleaned or easy autoclave for more then 4000 times. This contributes mostly to the elimination of cross-infection. VETSKOP® Series delivers the best of it with low friction with smooth-line chrome plating finish and ultra fine smooth edges. VETSKOP® Series bring a vast range of blades such as;

- RABBIT BLADE
- LEFT HAND MACINTOSH
- LEFT HAND MILLER WITH MAGNIFICATION
- LEFT HAND FLEXI-TIP
- LEFT HAND MILLER
- LEFT HAND MILLER (BAYONET) WITH MAGNIFICATION

AUG VETSKOP® Series are build with prime quality German made fiber. Each blade has more then 8000 strands with dia of 4.5mm, which delivers matchless illumination and bright shine intensity. The fiber is supremely polished from both ends to make if more clear.

### RABBIT BLADE

Our Rabbit Blade are specially designed to aid and simplify airway management in exotics, companion animals.

### LEFT HAND MACINTOSH

Our Left Hand Macintosh Blades are simply designed for daily basis clinical diagnose which help to elevate the epiglottis and to expose the laryngeal inlet, effortlessly.

### LEFT HAND MILLER

Our Left Hand Miller Blades are designed with smooth straight flange which helps to obtain a perfect view of the vocal cords by directly lifting the epiglottis or to expose the laryngeal inlet, effortlessly.

### LEFT HAND MILLER WITH MAGNIFICATION

Our Left Hand Miller Blades are designed with smooth straight flange which helps to obtain a perfect view of the vocal cords by directly lifting the epiglottis or to expose the laryngeal inlet, effortlessly. 3x Magnifying Lens with two-way swivel gives special focal distance.

### LEFT HAND FLEXI-TIP

Our Left Hand FLEXI-TIP Blades are designed with finest quality tilts tips with adjustment angle of 90° by pushing out the lever at back, the blades are most demanded for difficult intubation.

### LEFT HAND MILLER (BAYONET) WITH MAGNIFICATION

Our Left Hand Miller Blades are designed with smooth straight flange which helps to obtain a perfect view of the vocal cords by directly lifting the epiglottis or to expose the laryngeal inlet, effortlessly. Bayonet fitting. 3x Magnifying Lens.





## VETSKOP® | F.O RABBIT BLADE

### Features:

- Enclosed fiber optic core with 3.00mm diameter.
- 3000 plus fiber strands.
- High quality fiber ensuring brilliant light transmission.
- No sharp edges. Low profile straight blades.
- Specially designed for veterinarian.
- Free maintenance.
- Autoclave able upto 134 °C approximately 4000 times.
- No adverse effect in using ETO, steam autoclaving or by soaking in Sparax or Cidex.
- In accordance with ISO 7376/3.
- Elegant design with low profile heel with chrome plated finish.

REFERENCES	SIZE
------------	------

60.900.101	1
------------	---

## VETSKOP® | F.O LEFT HAND MACINTOSH BLADE

### Features:

- Enclosed fiber optic core with 4.5mm diameter.
- 8000 plus fiber strands.
- High quality fiber ensuring brilliant light transmission.
- No sharp edges. Low profile curved blades.
- Specially designed for veterinarian.
- Free maintenance.
- Autoclave able upto 134 °C approximately 4000 times.
- No adverse effect in using ETO, steam autoclaving or by soaking in Sparax or Cidex.
- In accordance with ISO 7376/3.
- Elegant design with low profile heel and chrome plated finish.



REFERENCES	SIZE	REFERENCES	SIZE
------------	------	------------	------

60.900.200	0	60.900.203	3
60.900.201	1	60.900.204	4
60.900.202	2	60.900.205	5



**AUG**  
MEDICAL

Quality. Defined



## VETSKOP® | F.O LEFT HAND FLEXI-TIP BLADE

### Features:

- Enclosed fiber optic core with 4.5mm diameter.
- 8000 plus fiber strands.
- High quality fiber ensuring brilliant light transmission.
- No sharp edges. Low profile curved blades.
- Specially designed for veterinarian.
- Free maintenance.
- Autoclave able upto 134 °C approximately 4000 times.
- No adverse effect in using ETO, steam autoclaving or by soaking in Sparax or Cidex.
- In accordance with ISO 7376/3.
- Elegant design with low profile heel with chrome plated finish.

REFERENCES	SIZE	REFERENCES	SIZE
60.900.401	1	60.900.403	3
60.900.402	2	60.900.404	4

## VETSKOP® | F.O LEFT HAND MILLER BLADE

### Features:

- Enclosed fiber optic core with 3.00mm diameter.
- 3000 plus fiber strands.
- High quality fiber ensuring brilliant light transmission.
- No sharp edges. Low profile straight blades.
- Specially designed for veterinarian.
- Free maintenance.
- Autoclave able upto 134 °C approximately 4000 times.
- No adverse effect in using ETO, steam autoclaving or by soaking in Sparax or Cidex.
- In accordance with ISO 7376/3.
- Elegant design with low profile heel with chrome plated finish.

REFERENCES	SIZE	REFERENCES	SIZE
60.900.299	00	60.900.302	2
60.900.300	0	60.900.303	3
60.900.301	1	60.900.304	4
		60.900.305	12"





## VETSKOP® | F.O LEFT HAND MILLER BLADE WITH MAGNIFICATION

### Features:

- Enclosed fiber optic core with 3.00mm diameter.
- 3000 plus fiber strands.
- High quality fiber ensuring brilliant light transmission.
- No sharp edges. Low profile straight blades.
- Specially designed for veterinarian.
- Free maintenance.
- Autoclave able upto 134 °C approximately 4000 times.
- No adverse effect in using ETO, steam autoclaving or by soaking in Sparax or Cidex.
- In accordance with ISO 7376/3.
- Elegant design with low profile heel with chrome plated finish.
- 3x Magnifying lens with two-way swivel gives special focal distance.

REFERENCES	SIZE	REFERENCES	SIZE
60.900.499	00	60.900.502	2
60.900.500	0	60.900.503	3
60.900.501	1	60.900.504	4
		60.900.505	12"

## VETSKOP® | F.O LEFT HAND MILLER BLADE WITH MAGNIFICATION (BAYONET FITTING)

### Features:

- High quality 3.00mm fiber ensuring brilliant light transmission.
- 2.5V LED is constructed with advanced German technology.
- Polished metal cap around lamp gives maximum supply of illumination.
- Lamp can be easily supersede with help of fingers.
- 50,000+ hours life of LED lamp with more than 130,000 Lux.
- No sharp edges. Low profile straight blades, Free maintenance.
- Autoclave able upto 134 °C approximately 4000 times.
- No adverse effect in using ETO, steam autoclaving or by soaking in Sparax or Cidex.
- Bayonet fitting.
- In accordance with ISO 7376/3.
- Elegant design with low profile heel with chrome plated finish.
- 3x Magnifying lens with two-way swivel gives special focal distance.



REFERENCES	SIZE	REFERENCES	SIZE
60.900.599	00	60.900.602	2
60.900.600	0	60.900.603	3
60.900.601	1	60.900.604	4
		60.900.605	12"





## VETSKOP<sup>®</sup> | F.O LED/XENON HANDLE

### Features:

- 2.5V LED/Xenon is constructed with advance German technology.
- Polished metal cap around Lamp gives maximum supply of illumination.
- 50+ hours life of LED/Xenon Lamp.
- More than 11,000 Lux.
- Closed battery cover prevent liquids penetration inside the battery compartment.
- Can be easily operated with 2 x C size dry cell batteries.
- Green color coding indicate universal fitting & standard of ISO 7376/3.
- Build in Brass.Premium chrome plating finishing.

REF	SIZE	LAMP	LENGHT X WIDTH
90.100.201	SMALL	LED	154mm x 19mm
90.100.202	MEDIUM	LED	153mm x 29mm
90.100.401	SMALL	LED RC	154mm x 19mm
90.100.402	MEDIUM	LED RC	153mm x 29mm
90.100.301	SMALL	XENON	154mm x 19mm
90.100.302	MEDIUM	XENON	153mm x 29mm
90.100.411	SMALL	XENON RC	154mm x 19mm
90.100.412	MEDIUM	XENON RC	153mm x 29mm


## VETSKOP<sup>®</sup> | POWER SOURCE RHEOSTAT HANDLE (BAYONET)

### Features:

- Light strength can be easily regulated with dimmable switch.
- Polished metal bayonet fitting gives maximum supply for illumination.
- Knurled/Grooved handle design provides excellent grip
- Closed battery cover prevent liquids penetration inside the battery compartment.
- Can be easily operated with 2 x C size dry cell batteries.
- Bayonet fitting based..
- Build in brass with premium chrome plating finishing.

REF	SIZE	POWER	LENGHT X WIDTH
60.700.941	SMALL	AA	133mm x 18mm
60.700.951	SMALL	AA (RC)	133mm x 18mm
60.700.940	MEDIUM	C	143mm x 29mm
60.700.950	MEDIUM	C (RC)	149mm x 29mm






## VETSKOP<sup>®</sup> | OTOSCOPE

Elegant Veterinary Otoscope Designed With An All Metal Housing & Scratch Resistant Glass Lenses Which Guarantee A Lifetime Of Maintenance-free Use And Improve The Efficiency Of Examinations & Diagnoses. Get A Wider, Deeper, Clearer View Of The Ear With The AUG<sup>®</sup> VETSKOP<sup>®</sup>, An Advance In Hand-held Otoscope That Helps You Perform One Of Your Most Common Procedures, Thank To The Unique German LED Technology Gives Concentrated Bright Light For Perfect Illumination. Dimmable Switch Helps You To Regulate The Light. High Quality Metal Handle With Bayonet Fitting & Knurled Surface Makes The Product Grip Secured And Efficient. AUG<sup>®</sup> VETSKOP<sup>®</sup> Otoscope Has Excellent Quality And High Durability With Long Lasting Companion For Veterinary Diagnostics.

## VETSKOP<sup>®</sup> | OPHTHALMOSCOPE

Elegant Veterinary Ophthalmoscope Designed With A Black Plastic Head Housing And Excellent Optics With Best Color Rendering In Minimum Dimensions, Ideal For Rounds And House Calls. AUG<sup>®</sup> Special Optics Produce Crisp, Clear Images Of The Whole Area Under Examination And Minimize Reflections Which Guarantee A Lifetime Of Maintenance-free Use And Improve The Efficiency Of Examinations & Diagnoses. High Quality Bayonet Fitting Metal Handle With Knurled Surface Which Makes The Product Grip Secured And Efficient. AUG<sup>®</sup> VETSKOP<sup>®</sup> Ophthalmoscope Has Excellent Quality And High Durability With Long Lasting Companion For Comprehensive Pricing. Unique German XHL Vacuum Technology Gives Concentrated Bright Light For Perfect Illumination With An Advance Feature Of Dimmable Switch.





**REFERENCES (In Case)**

20.000.600

**VETSKOP® | F.O LED OTOSCOPE**

**FEATURES**

- Enclosed fiber optic core with 4.00mm diameter
- High quality fibers ensuring brilliant light transmission.
- 2.5V German LED technology gives brighter illumination.
- 3x Magnifying Lens with two-way swivel gives special focal distance.
- Light strength can be easily regulated with dimmable switch.
- Operated with 2 x C size batteries.
- Premium chrome plated finishing.
- Reusable Metal Ear Funnels Dia Ø 2.5mm, 3.5mm, 4.5mm, 5.5mm & 9.4mm
- Compatible Vaginal Speculum is also available.
- In Plastic Case with Spare Lamp.
- Easy replacement of bulb.



**REFERENCES (In Case) WSB**

20.000.605



**REFERENCES (In Case)**

20.000.601

**VETSKOP® | F.O XENON OTOSCOPE**

**FEATURES**

- Enclosed fiber optic core with 4.00mm diameter
- High quality fibers ensuring brilliant light transmission.
- 2.5V German Xenon technology gives brighter illumination.
- 3x Magnifying Lens with two-way swivel gives special focal distance.
- Light strength can be easily regulated with dimmable switch.
- Operated with 2 x C size batteries.
- Premium chrome plated finishing.
- Reusable Metal Ear Funnels Dia Ø 2.5mm, 3.5mm, 4.5mm, 5.5mm & 9.4mm
- Compatible Vaginal Speculum is also available.
- In Plastic Case with Spare Lamp.
- Easy replacement of bulb.



**REFERENCES (In Case) WSB**

20.000.606





## VETSKOP® | VACUUM OTOSCOPE

### FEATURES

- 2.5V German XHL Vacuum technology gives brighter illumination.
- 3x Magnifying Lens with two-way swivel gives special focal distance.
- Light strength can be easily regulated with dimmable switch.
- Operated with 2 x C size batteries.
- Premium chrome plated finishing.
- Reusable Metal Ear Funnels Dia Ø 2.5mm, 3.5mm, 4.5mm, 5.5mm & 9.4mm
- Compatible Vaginal Speculum is also available.
- In Plastic Case with Spare Lamp.
- Easy replacement of bulb.



### REFERENCES (In Case)

20.000.602

### REFERENCES (In Case) WSB

20.000.607



## VETSKOP® | LED OPHTHALMOSCOPE

### FEATURES

- 2.5V German LED Lamp gives brighter illumination.
- 5 Different Apertures with 2 RED-FREE filters & 1 neutral options;
 

⏏	○	⊛	○	◐	&	●	●	●
---	---	---	---	---	---	---	---	---
- **Small Spot:** For reduced reflection with small pupils.
- **Large Spot:** For normal funduscopy.
- **Semi-Circle:** For reduced reflection with small pupils.
- **Fixation Star:** For determination of central or eccentric fixation, suitable for examining children.
- **Green or Blue Filter:** For contrast enhancement when assessing six vascular disorders.
- **Slit or Streak:** Helpful in determining various levels of lesions, particularly tumors and edematous disc.
- Color coded lens for excellent resolution:
 

+	0	1	2	3	4	5	6	7	8	9	10	12	15	20	40	
-		1	2	3	4	5	6	7	8	9	10	15	20	25	30	35
- Light strength can be easily regulated with dimmable switch.
- Operated with 2 x C size batteries.
- Bayonet fitting based.
- Premium chrome plated finishing.
- In Fabric Pouch.
- Also available in hard plastic case with spare lamp.
- Easy replacement of bulb.



### REFERENCES (In Pouch)

20.000.620



### REFERENCES (In Case) WSB

20.000.625





## VETSKOP® | XENON OPHTHALMOSCOPE

### FEATURES

- 2.5V German Xenon Lamp gives brighter illumination.
- 5 Different Apertures with 2 RED-FREE filters & 1 neutral options;
 

	○	⊕	○	◐	&	●	●	●
--	---	---	---	---	---	---	---	---
- **Small Spot:** For reduced reflection with small pupils.
- **Large Spot:** For normal funduscopy.
- **Semi-Circle:** For reduced reflection with small pupils.
- **Fixation Star:** For determination of central or eccentric fixation, suitable for examining children.
- **Green or Blue Filter:** For contrast enhancement when assessing six vascular disorders.
- **Slit or Streak:** Helpful in determining various levels of lesions, particularly tumors and edematous disc.
- Color coded lens for excellent resolution:
 

+	0	1	2	3	4	5	6	7	8	9	10	12	15	20	40	
-		1	2	3	4	5	6	7	8	9	10	15	20	25	30	35
- Light strength can be easily regulated with dimmable switch.
- Operated with 2 x C size batteries.
- Bayonet fitting based.
- Premium chrome plated finishing.
- In Fabric Pouch.
- Also available in hard plastic case with spare lamp.
- Easy replacement of bulb.



### REFERENCES (In Pouch)

20.000.621 

### REFERENCES (In Case)

20.000.626 

## VETSKOP® | RC LED/XENON OPHTHALMOSCOPE

### FEATURES



- 3.5V German LED/Xenon Lamp gives brighter illumination.
- 5 Different Apertures with 2 RED-FREE filters & 1 neutral options;
 

	○	⊕	○	◐	&	●	●	●
--	---	---	---	---	---	---	---	---
- **Small Spot:** For reduced reflection with small pupils.
- **Large Spot:** For normal funduscopy.
- **Semi-Circle:** For reduced reflection with small pupils.
- **Fixation Star:** For determination of central or eccentric fixation, suitable for examining children.
- **Green or Blue Filter:** For contrast enhancement when assessing six vascular disorders.
- **Slit or Streak:** Helpful in determining various levels of lesions, particularly tumors and edematous disc.
- Color coded lens for excellent resolution:
 



+	0	1	2	3	4	5	6	7	8	9	10	12	15	20	40	
-		1	2	3	4	5	6	7	8	9	10	15	20	25	30	35
- Light strength can be easily regulated with dimmable switch.
- Operated with AUG® 3.5V Ni-MH 1000mAh Rechargeable Battery (Included).
- Can be easily recharge with AUG® Pro Charger.
- 2 hour charge gives more than 3 hours working capacity with 3.5V Ni-MH 1000mAh rechargeable battery.
- Premium chrome plated finishing.
- Available in pouch, also in hard plastic case with spare lamp.
- Easy replacement of bulb.



### REFERENCES (In Pouch)

20.000.630   
 20.000.631 

### REFERENCES (In Case) WSB




20.000.635   
 20.000.636 





### VETSKOP® | OTO-OPHTHALMOSCOPE SET

#### REFERENCES (In Case) WSB


20.000.650		20.000.656	
20.000.651		20.000.657	
20.000.652		20.000.658	

\*\*See feature and details on page 30-31



### VETSKOP® | OTO-OPHTHALMOSCOPE SET

#### REFERENCES (In Case) WSB







20.000.670		20.000.676	
20.000.671		20.000.677	
20.000.672		20.000.678	

\*\*See feature and details on page 30-31



### VETSKOP® | DIAGNOSTIC SET

#### REFERENCES (In Case) WSB

20.000.680		20.000.686	
20.000.681		20.000.687	
20.000.682		20.000.688	

\*\*See feature and details on page 30-31



## DESKTOP CHARGER

### RECHARGE Ni-MH & Li-ON BATTERIES

#### FEATURES

- Charge two handles simultaneously.
- Arc shaped charging indicator with breathing effects.
- Compatible of Ni-MH & Li-ON batteries.
- Compatible with all rechargeable handles of PRO Series.
- Input 220V, Takes 2.5 hours to recharge 1000 mAh battery.



## USB 2.0 CHARGER

### RECHARGE Li-ON BATTERIES 5V-1.0A

#### FEATURES

- Charge handle simultaneously while instrument in operation.
- Green / Red charging indicator for battery power.
- Takes 2 hours to recharge 1000 mAh battery.
- Compatible with all rechargeable handles of PRO Series.



## BATTERY PACK

### RECHARGEABLE Li-ON & Ni-MH

#### FEATURES

- Support all PRO Series Handles.
- Easy to replace with protective packing.
- Takes 2 hours to recharge 1000 mAh battery.
- Compatible chargers are Desktop & USB2.0

VOLT	REFERENCES	DETAIL
2.5V/3.5V	<a href="#">60.300.100</a>	DESKTOP CHARGER
2.5V/3.5V	<a href="#">60.300.150</a>	USB 2.0 CHARGER
2.5V	<a href="#">90.100.110</a>	Ni-MH SMALL BATTERY
2.5V	<a href="#">90.100.115</a>	Li-ON SMALL BATTERY
3.5V	<a href="#">90.100.120</a>	Ni-MH MEDIUM BATTERY
3.5V	<a href="#">90.100.125</a>	Li-ON MEDIUM BATTERY

# ACCESSORIES



## REUSABLE EARTIP (Metal)

△ Closed Metal tip	Ref:	△ Open Gutter Metal tip	Ref:
57.0 x 4.0mm	60.200.701	65.0 x 6.0mm	60.200.704
65.0 x 7.5mm	60.200.702	90.0 x 7.0mm	60.200.705
90.0 x 9.0mm	60.200.703		



## REUSABLE VAGINAL SPECULUM (Metal)

### REFERENCES

60.200.750

## LAMPS FOR LARYNGOSCOPE



LED LAMP

REF: 60.100.300 Volt: 2.5V

REF: 60.100.310 Volt: 3.5V



XENON LAMP

REF: 60.100.450 Volt: 2.5V

REF: 60.100.451 Volt: 3.5V



VACUUM LAMP

REF: 60.100.400 SIZE: SMALL

REF: 60.100.410 SIZE: LARGE

## LAMPS FOR OTOSCOPE



GERMAN LED LAMP

REF: 60.100.200 Volt: 2.5V

REF: 60.100.205 Volt: 3.5V



GERMAN XENON LAMP

REF: 60.100.600 Volt: 2.5V

REF: 60.100.605 Volt: 3.5V



GERMAN VACUUM LAMP

REF: 60.100.650 Volt: 2.5V

REF: 60.100.655 Volt: 3.5V

## LAMPS FOR OPHTHALMOSCOPE



GERMAN LAMP | LED & XENON

REF: 60.100.700 Light: LED

REF: 60.100.750 Light: Xenon



GERMAN LAMP | LED & XENON

REF: 60.100.915 Light: LED 2.5V

REF: 60.100.905 Light: Xenon 2.5V

REF: 60.100.910 Light: LED 3.5V

REF: 60.100.900 Light: Xenon 3.5V



GERMAN VACUUM LAMP

REF: 60.100.090 Volt: 2.5V

REF: 60.100.091 Volt: 3.5V



**HEADQUARTER:**  
**AUG Medical Sp. z o.o.**  
Ul. Poleczki 21 Ramp 13  
02-822 Warszawa-Poland  
Tel: +48 22 243 37 09  
Email: [info@augmedical.pl](mailto:info@augmedical.pl)  
Email: [info@augmedical.com](mailto:info@augmedical.com)

**AUG MEDICAL LLC**  
777 Brickell Ave. Suite 500  
Miami Florida 33131 USA  
Tel: +1 305 671 9396  
Email: [info@augmedical.com](mailto:info@augmedical.com)

